SONOPULS 492 NEW MODEL



The top device of the 4 series is the Sonopuls 492. This complete unit ensures that you can apply ultrasound therapy, electrotherapy and combination therapy. You will have a solution for every disorder at your fingertips with both 1 MHz and 3 MHz ultrasound, with a small and/or a large ultrasound treatment head and with all modern and effective current types that are fully integrated!

The neatly arranged touch screen operating panel makes operation extremely easy and fast. All information that you may need is displayed in a very well-organised manner. Selecting the correct parameters is also easy. Thanks to the integrated evidence-based protocols, the best recommendations will be made. The graphic help screens that display localisations and anatomic representations guarantee the correct execution of the therapy.

The Sonopuls 492 can be used as a mobile device thanks to the battery (optional). In addition to using it with a mains power supply, the equipment can also be used in a clinic, next to the bed of a patient or at home. Thus the multifunctional Sonopuls 492 can also become your assistant wherever you may go.

The Sonopuls 492 can be expanded with the Vacotron 460 suction module. With this Vacotron, that can be integrated beautifully with the Sonopuls 492, you can easily position electrodes. A combination that will make your job even easier and faster to accomplish.

The Sonopuls 492: the complete, easy-to-use, fast and responsible choice for ultrasound therapy, electrotherapy and combination therapy!



ULTRASOUND THERAPY

Equipment for therapy using high frequency sound waves. The high frequency vibrations can cause a mechanical and thermal effect in the tissue. Ultrasound therapy promotes the body's capacity for selfhealing and tissue repair.

ELECTROTHERAPY

Equipment for pain management, muscle stimulation and/or diagnostics. Electrotherapy applications for fast pain relief without side effects. Muscles can be stimulated to improve the function of soft tissue. Electrotherapy can be used to measure differentiated tissue response to stimulation, to provide a diagnosis about tissue quality.

COMBINATION THERAPY

Equipment for electrotherapy and ultrasound therapy. This type of therapy can be applied side-by-side and at the same time, providing a very efficient therapy. Pain reduction in combination with ultrasound therapy provides rapid post-traumatic treatment, resulting in tissue repair in a minimum of time.

SONOPULS 492

1498913 SONOPULS 492 WITH 2 MULTI-FREQUENCY TREATMENT HEADS (5 CM² AND 0,8 CM²)



1498911 SONOPULS 492 WITH MULTI-FREQUENCY TREATMENT HEAD, LARGE (5 CM²)



1498912 SONOPULS 492 WITH MULTI-FREQUENCY TREATMENT HEAD, SMALL (0,8 CM²)



NEW MODEL

TECHNICAL SPECIFICATIONS

Power supply Dimensions (WxDxH) Device Device with base Device with Vacotron Weight Weight, including the optional battery $100 - 240 \text{ V} \pm 10\% \text{ / } 50 \text{ / } 60 \text{ Hz}$

24 x 32 x 12 cm

2 kg

3 kg

2

5 cm²

0,8 cm²

1000 +

1000+

1 en 3 MHz

24 x 30.5 x 18.2 cm

24 x 30.5 x 21.6 cm

continuous and pulsed 5, 10, 20, 33, 50, 80%

16 Hz, 48 Hz en 100 Hz

0-2 W/cm² continuous,

0-3 W/cm² pulsed⁽¹⁾

25 - Evidence Based

2 (100% independent)

42 - Evidence Based

36 (per channel)

ULTRASOUND

Ultrasound frequencies Ultrasound Duty cycles Pulse frequency Intensity

Number of US connections ERA US applicator large ERA US applicator small Pre-programmed treatment suggestions Programmable positions

ELECTROTHERAPY

Current channels Current types Pre-programmed treatment suggestions Programmable positions

COMBINATIONTHERAPY

Pre-programmed treatment suggestions 2 - Evidence Based Programmable positions 1000+

⁽¹⁾ = Except 80%, which has a maximum of 2,5 W/cm²

14987574-Series Operating Instructions (cd rom)14987564-Series Information Booklet3444290Mains cable 230V-EUR1498010Device base (inclination support) (not for Vacotron 460)xxxxxxUltrasound multi-frequency treatment head(s) (model dependent)xxxxxxxHolder treatment head (model dependent)3442929(1)Contact-gel, bottle 250 ml, 1x3444021Strap 250x3 cm3444020Strap 100x3 cm

- 2x 3444129 Rubber electrodes 6x8 cm, 2 mm female, set of 2 pcs.
- 1460266Moist pads for rubber electrodes 6x8 cm, set of 4 pcs.2x3444211Patient cable 2-core & 2 mm male plugs black,
 - with coloured clips

STANDARD ACCESSORIES SONOPULS 492

 $^{(1)}$ = The Sonopuls is delivered with 1 bottle of contact-gel. The articlenumber 3442929 however represents a box of 12 bottles.

> ENRAF-NONIUS P.O. Box 12080 3004 GB Rotterdam The Netherlands T.: +31-(0)10 - 203 06 00 E: info@enraf-nonius.nl www.enraf-nonius.com

SONOPULS 492

NEW MODEL

EVIDENCE-BASED PROTOCOLS

- · 25 Evidence-Based treatment protocols for ultrasound,
- 42 Evidence-Based protocols for electrotherapy,
- · 2 Evidence-Based protocols for combination therapy

Complete including the display of the anatomic localisations and the treatment instructions

Clinical Protocols Clinic



The patient can be involved in his or her own treatment by using the optional "remote control". Especially when using electrotherapy, the patient will have the possibility to control the therapy

REMOTE CONTROL (OPTIONAL)

ERGONOMIC DESIGN

The treatment head is made to fit your hand comfortably. A design completely based on functionality and user friendliness.



BATTERY-OPERATED (OPTIONAL) AND, THERE-FORE, MOBILE

In addition to being able to use the 4-series unit connected to a mains power supply, the device can also be used with a batteryoperated power supply. This will ensure that you can go everywhere with your unit, regardless of whether this is at the patient's home or at his or her hospital bed.





SONOPULS 492



MADE BY ENRAF-NONIUS

Made by Enraf-Nonius and, therefore, safe, reliable and durable. A device made by Enraf-Nonius, THE supplier for the physiotherapist worldwide. Where quality, reliability and durability come as standard.

NEW MODEL

CRYSTAL CLEAR TFT COLOUR SCREEN

A crystal clear, full colour and scratch resistant screen that allows full visibility on all parameters from all viewing angles.



TOUCH SCREEN USER INTERFACE

Easy to use by means of a clear menu via a touch screen interface. Extremely intuitive and as easy as a state-of-the-art smart phone.

COMPACT AND STYLISH

An extremely compact device, small and light-weight for easy transportation. At the same time stylish and beautiful, a device that is worth-while showing!



2 TREATMENT HEADS POSSIBLE

You can connect both the small 0,8 cm2 as well as the large 5 cm² ultrasound treatment head in order to quickly change heads.

SONOPULS 492

NEW MODEL

PROGRAMMED FOR SPEED

The menu structure is optimized to speed. In most cases you can start a pathology-specific treatment in just 3 steps! Easy, quick and simple.





CONTACT CONTROL

A 360 degrees light ring for contact control. A clear indication of the effective transfer of energy that enables you to conduct the therapy as efficient as possible.

A TREATMENT HEAD WITH ADDED VALUE

An intelligent and ergonomic treatment head, packed with electronics and covered with a special coating making it smooth and comfortable to use.



| P Electrotherapy | Favorite name |
|-----------------------|----------------------------|
| | MR ATKINSON TENS 2 CHANNEL |
| Clinical Protocols | 1 2 3 4 6 6 7 8 9 0 |
| Favorites | qwertyulop asdfohjki |
| | 🛆 Z X C V b n m 🔒 - |
| Manual Operation | All¢ |
| Programming | |
| | |
| | |
| ■ 0:00 🦞 ■ 0:00 🐔 ■ 5 | |
| B 0:00 👎 B 0:00 🛃 S | ♥ ■ 0.00 ♥ ■ ≦^ |

PROGRAMMABLE


In addition to the standard treatment protocols, you can also save your personal protocols (including sequential applications) based on your own guidelines.

ACCESSORIES SONOPULS 492

| JLTRASO | UND ACCESSORIES | | <u>a</u> | | |
|---|---|--------------------|-------------------------|----------------------------------|----------------|
| | UND CONTACT-GEL | | | | S |
| 3442929 | Contact-gel, bottle 250 ml, box of 12 | | | | |
| 3442930 | Contact-gel, bottle 850 ml, box of 12 | 3442929 | 3442930 | 3442931 | 3442932 |
| 3442931 | Contact-gel, canister of 5 I | | | | |
| 3442932 | Dispenser-set for 5 I canister | | | | |
| | | | | | 5 |
| | DES AND MOIST PADS E ELECTRODES | | ** | 00 | |
| 3444222 | Adhesive electrodes Ø 2,2 cm, 2 mm female, 10 sheets of 8 (also for EMG) | 3444222 | 3444056 | 3444135 | 3444057 |
| 3444056 | EN-Trode Ø 3,2 cm, 2 mm female,10 sheets of 4 | | | | |
| 3444135 | EN-Trode Ø 5,0 cm, 2 mm female, 10 sheets of 4 | | | | |
| 3444057 | EN-Trode 5x5 cm, 2 mm female, 10 sheets of 4 | | | | |
| 3444058 | EN-Trode 5x9 cm, 2 mm female, 10 sheets of 4 | | | | |
| | ELECTRODES | 3444058 | | | |
| 3444128 3444129 | Rubber electrodes 4x6 cm, 2 mm female, set of 2. Rubber electrodes 6x8 cm, 2 mm female, set of 2 | | | | |
| 3444129 | Rubber electrodes 6x8 cm, 2 mm female, set of 2 Rubber electrodes 8x12 cm, 2 mm female, set of 2 | 3 | | | |
| | DS FOR RUBBER ELECTRODES | | | | |
| 460273 | Moist pads for rubber electrodes 4x6 cm, set of 4 | 3444128 | 3444129 | 3444130 | |
| 460266 | Moist pads for rubber electrodes 6x8 cm, set of 4 | | | 00177700 | |
| 460275 | Moist pads for rubber electrodes 8x12 cm, set of 4 | | | | |
| IXATION | · | | | | |
| 3444020 | Strap 100x3 cm | | | | |
| 3444021 | Strap 250x3 cm | 1460273 | 1460266 | 1460275 | |
| /////// | Strap 100x5 cm | | | | |
| 3444022 | | | | | - |
| | Strap 250x5 cm | | | | |
| 3444022 3444023 | • | | 00 | 00 | 00 |
| 3444022 3444023 | Strap 250x5 cm | | | 2 3 4 | 3444023 |
| 3444022 3444023 POINT EL | Strap 250x5 cm ECTRODES | € 3444020 | 3444021 | 3444022 | ¢ 3444023 |
| 3444022 3444023 POINT EL | Strap 250x5 cm ECTRODES Electrode holder for 1480801 Point electrode 5 mm, 4 mm female (in combination with | 3444020 | 3444021 | 3444022 | 4 44023 |
| 3444022 3444023 POINT EL 1480800 1480801 1460369 | Strap 250x5 cm ECTRODES Electrode holder for 1480801 Point electrode 5 mm, 4 mm female (in combination with adapters 2523523 / 2523524) Point electrode, 17 mm, 4 mm female (in combination with adapters 2523523 / 2523524) | 3444020 1480800 | € 3444021 1480801 | 3444022 3444022 1460369 | 3444023 |
| 3444022 3444023 POINT EL 1480800 1480801 1460369 | Strap 250x5 cm ECTRODES Electrode holder for 1480801 Point electrode 5 mm, 4 mm female (in combination with adapters 2523523 / 2523524) Point electrode, 17 mm, 4 mm female (in combination with adapters 2523523 / 2523524) S | | ~ | . | 3444023 |
| 3444022 3444023 POINT EL 1480800 1480801 1460369 | Strap 250x5 cm ECTRODES Electrode holder for 1480801 Point electrode 5 mm, 4 mm female (in combination with adapters 2523523 / 2523524) Point electrode, 17 mm, 4 mm female (in combination with adapters 2523523 / 2523524) | | ~ | . | € 3444023 |

01-2014- Enraf-Nonius Electr





OTHER PRODUCTS THAT MIGHT INTEREST YOU

VACOTRON 460

The Sonopuls 492 can be extended with a vacuum module, the Vacotron 460. A vacuum unit for easy application of vacuum electrodes for combination therapy and electrotherapy.

EN-CAR U3

The EN-Car U3 is your companion for any piece of Enraf-Nonius equipment. For laser, ultrasound, biofeedback, shockwave or electrotherapy, with the EN-Car U3 you have it all close at hand.

ORDERING DATA

1498950 Vacotron 460 1468960 EN-Car U3

Also see the seperate product leaflets at WWW.ENRAF-NONIUS.COM/ELECTROTHERAPY





